



Biometrics  
Quality of life

MOLECULAR AND  
CLINICAL DIAGNOSTICS

# NucleoPURE

NUCLEIC ACID EXTRACTION SYSTEM



[www.biometrics-technologies.com](http://www.biometrics-technologies.com)

# NucleoPURE Extraction System

## Automated Nucleic Acid Extraction System

With the latest isolation method of nucleic acid transferred by magnetic rods and preloaded magnetic beads-based extraction reagent, The NucleoPURE Nucleic Acid Automatic Extraction System can process up to 96 samples and automatically isolate nucleic acids from various samples such as blood, cells, viruses, etc. The magnetic beads are adsorbed, transferred and released by the special magnetic rod, so as to realize the fully automated purification of nucleic acid.



### Advanced Stirring Technology

Stirring speed up to 10HZ (10 times up and down in 1 second) enables effective mixing of samples with high viscosity (blood, feces, etc.).




### Safe and Reliable

The operation will automatically suspend once the safety door opens, using a front-door sensor system, and will only resume after the safety door is shut, providing enhanced protection for users.



### Self-Diagnostics

Auto checks before every run ensures that the sleeve, rod and well plate are positioned correctly thereby avoiding any errors due to omission or interference.



### Multi Heating Modules

Multiple heating modules with individual heating controls ensures temperature accuracy for lysis and elution.



### Stable and Efficient

Automated processing avoids discrepancies and errors caused by manual operations and increases the efficiency




### Ease of Use

A 10-inch color touch screen and user-friendly interface lets you easily set purification parameters.



### Contamination Control

Built-in UV disinfection lamp eliminates cross contamination



### Compact & Light Weight

NucleoPURE 32 help to save your limited laboratory bench space



# NucleoPURE 32/96 EXTRACTION SYSTEM


## BUILT-IN SOFTWARE

### Run Status Monitoring

The **Run Status** interface shows real-time extraction progress, active wells, temperatures, timers, with controls like **Start**, **Pause**, and **End** for user interaction.





**NucleoPURE 32**



**NucleoPURE 96**



### Personalized Shortcuts and Extensive Program Capacity

Enhance efficiency with the ability to create custom shortcuts directly on the screen for frequently used operations. Additionally, the system supports the creation and storage of over 500 extraction programs, providing unparalleled flexibility and convenience for diverse workflows.



### Smart Auto Check for Seamless Operation


With this automated check, you can confidently start your program, knowing everything is set up perfectly!





# NucleoPURE Extraction System

## Working Principle




- 1 Sample
- 2 The sample are lysed in lysis buffer to release the nucleic buffer
- 3 Transfer the magnetic beads into lysis buffer blend fully and make the nucleic acids absorbed onto the specific coated material on the surface of the magnetic beads
- 4 Transfer the magnetic beads to the elution buffer and mix thoroughly, the nucleic acid falls off the surface of the magnetic beads and dissolves into the elution buffer.
- 5 Clean the surface of magnetic beads repeatedly, to remove unnecessary protein, salt or other impurities

# NucleoPURE 32/96 EXTRACTION SYSTEM


## WORKING STEP

### SAMPLE LOADING




### LYSIS/NUCLEIC ACID BINDING

Lysis buffer solution promotes cell rupture and releases nucleic acids. Magnetic beads transferred from wells 2/8 to the lysis wells binds to the target nucleic acid molecules. The magnetic beads with the bound particles are then moved to the wells 3/9 for washing.




### WASHING

Alcoholic wash buffers in the wash wells facilitate the removal of contaminants. Since the magnetic sleeves were soaked in alcohol during the wash stage, they need to 'air dry' above the wells so that the alcohol can volatilize before starting the elution stage



### ELUTION


After ethanol is removed, the nucleic acid bound to the magnetic beads are mixed in elution buffer, releasing the nucleic acid.



## SPECIFICATION



<b>Model</b>	<b>NucleoPURE 32</b>
<b>Volume Capacity</b>	50-1000 $\mu$ L
<b>Sample Quantity</b>	1-32 samples
<b>Plate Quantity</b>	2 plates
<b>Plate Type</b>	96-well deep well plate
<b>Processing Time</b>	20-40 min/run (depending on reagent types and volumes)
<b>Reagent Type</b>	Magnetic bead-based open source reagent
<b>Program Management</b>	Create, Edit, and Delete options available
<b>Magnetic Force</b>	5000 gauss
<b>Heating Temp.</b>	Room temp. - 125°C
<b>Shaking Mixing</b>	Adjustable: 8 Levels
<b>Max. Stirring Speed</b>	10 Hz
<b>Magnetic Bead Recovery Rate</b>	$\geq$ 96%
<b>Inter-well Purification Variation</b>	CV <3%
<b>Working Noise</b>	<55 dB
<b>Exhaust</b>	Fan
<b>Data Storage Capacity</b>	> 500 programs
<b>Dimension</b>	349 (W) x 325 (D) x 390 (H) mm
<b>Display</b>	10-inch color touch screen
<b>Lighting</b>	Yes
<b>Contamination Control</b>	UV light sterilization, Built-in HEPA filters (Option)
<b>Power Supply</b>	200-240 VAC, 100-120 VAC; 50/60 Hz
<b>Weight</b>	25 Kg
<b>Extension Ports</b>	2 USB port, 1 Ethernet port



## SPECIFICATION

<b>Model</b>	<b>NucleoPURE 96</b>
<b>Volume Capacity</b>	50-1000 $\mu$ L
<b>Sample Quantity</b>	1-96 samples
<b>Plate Quantity</b>	6 plates
<b>Plate Type</b>	96-well deep well plate
<b>Processing Time</b>	20-40 min/run (depending on reagent types and volumes)
<b>Reagent Type</b>	Magnetic bead-based open source reagent
<b>Program management</b>	Create, Edit, and Delete options available
<b>Magnetic Force</b>	5000 gaussses
<b>Heating Temp.</b>	Room Temp. - 125°C
<b>Shaking Mixing</b>	Adjustable: 8 Levels
<b>Max. Stirring Speed</b>	10 Hz
<b>Magnetic Bead Recovery Rate</b>	$\geq$ 96%
<b>Inter-well Purification Variation</b>	CV <5%
<b>Working Noise</b>	<55 dB
<b>Exhaust</b>	Fan
<b>Data Storage Capacity</b>	> 500 programs
<b>Dimension</b>	349 (W) x 460 (D) x 390 (H) mm
<b>Display</b>	10-inch color touch screen
<b>Lighting</b>	Yes
<b>Contamination Control</b>	UV light sterilization, Built-in HEPA filters (Option)
<b>Power Supply</b>	200-240 VAC, 100-120 VAC; 50/60 Hz
<b>Weight</b>	35 Kg
<b>Extension Ports</b>	2 USB port, 1 Ethernet port



# Biometrics

Quality of life

**Biometrics Technologies, Inc. (Headquarter)**

1220 N. Market St. Suite 806, Wilmington,  
DE 19801 USA

Email: [info@biometrics-technologies.com](mailto:info@biometrics-technologies.com)

Website: [www.biometrics-technologies.com](http://www.biometrics-technologies.com)

**Biometrics Technologies Co., Ltd. (Taiwan)**

15F-8, No. 99, Sec. 1, Xintai 5th Rd.,  
Xizhi District, New Taipei City 221, Taiwan.

Email: [sales@biometrics-technologies.com](mailto:sales@biometrics-technologies.com)

**Biometrics Technologies (Asia Pacific) Ltd.  
(International & Asia Pacific Support Center)**

18, 7<sup>th</sup> Fl. Sricharoenchai Bldg., Tiwanon Rd.,  
Talat Khwan, Mueang, Nonthaburi 11000 Thailand.

Email: [info@biometrics-technologies.com](mailto:info@biometrics-technologies.com)

