

## Determination of Total Polyphenol Content in Green Tea Extract by Folin–Ciocalteu Assay

### Microplate Reader Application | MR-1000S

#### ■ Application

Quantitative determination of total polyphenol content (TPC) in green tea extract using the Folin–Ciocalteu colorimetric assay. The assay is measured at 765 nm and quantified against a gallic acid calibration curve, with results expressed as gallic acid equivalents (GAE).

#### ■ Principle

Phenolic compounds react with Folin–Ciocalteu reagent in the presence of sodium carbonate to form a blue-colored complex. The color intensity increases with the amount of reactive phenolic compounds. Absorbance is measured at 765 nm and compared with gallic acid standards.

#### ■ Workflow

Green tea extract → Add Folin–Ciocalteu reagent → Add Na<sub>2</sub>CO<sub>3</sub> → Incubate in the dark → Read OD<sub>765</sub> → Gallic acid standard curve → Calculate TPC as GAE



#### ■ Instrument Setup

| Parameter        | Recommended setup                                                     |
|------------------|-----------------------------------------------------------------------|
| Instrument       | MR-1000S Microplate Reader                                            |
| Measurement mode | Endpoint absorbance                                                   |
| Wavelength       | 765 nm                                                                |
| Plate            | 96-well microplate                                                    |
| Replicates       | Standards, blank and samples in replicate                             |
| Quantitation     | Gallic acid calibration curve → unknown concentration → report as GAE |

#### ■ Material and Standards

| Item                    | Published method condition           |
|-------------------------|--------------------------------------|
| Green tea extract       | Prepared in water at 100 µg/mL       |
| Gallic acid standards   | 10, 20, 30, 40 and 50 µg/mL          |
| Blank                   | Water                                |
| Folin–Ciocalteu reagent | 10% solution                         |
| Sodium carbonate        | 7.5% Na <sub>2</sub> CO <sub>3</sub> |
| Item                    | Published method condition           |

#### ■ Microplate Assay Procedure

| Step                                        | Volume / condition          |
|---------------------------------------------|-----------------------------|
| 1. Add standard, sample or blank            | 50 µL per well              |
| 2. Add 10% Folin–Ciocalteu reagent          | 50 µL per well              |
| 3. Initial reaction                         | Stand for 3 min             |
| 4. Add 7.5% Na <sub>2</sub> CO <sub>3</sub> | 50 µL per well and mix      |
| 5. Incubation                               | 60 min at 25°C, in the dark |
| 6. Measurement                              | Read absorbance at 765 nm   |

## ■ Result

Based on the measured OD765 values, the microplate reader automatically calculated the concentrations of the unknown samples from the standard curve. Green tea B showed the highest calculated concentration (37.0 µg/mL), followed by Green tea A (25.2 µg/mL) and Beverage C (14.2 µg/mL), demonstrating quantitative determination of total polyphenol content using the gallic acid calibration method.




Figure 1: Gallic acid calibration curve

Table 1 Quantitative determination of total polyphenol content in unknown samples using the gallic acid calibration curve

| Sample      | OD765 | GAE Concentration |
|-------------|-------|-------------------|
| Green tea A | 0.402 | 25.2 µg/mL        |
| Green tea B | 0.565 | 37.0 µg/mL        |
| Beverage C  | 0.250 | 14.2 µg/mL        |

## ■ Conclusion

The MR-1000S enables quantitative total polyphenol analysis by measuring absorbance at 765 nm, generating a gallic acid calibration curve, and automatically calculating the concentrations of unknown samples. Results can be expressed as gallic acid equivalents (GAE), providing a convenient solution for routine food and beverage quality control.

## ■ Reference

AOAC Official Method 2017.13. Total Phenolic Content in Extracts: Folin-C (Folin and Ciocalteu) Colorimetric Method.  
 Comparison of colorimetric methods for the analysis of total polyphenols in green tea extracts. Bioscience, Biotechnology, and Biochemistry. 2024;88(7):798–805.  
 Determination of Total Phenolic Content Using the Folin-C Assay: Single-Laboratory Validation, First Action 2017.13. Journal of AOAC INTERNATIONAL. 2018;101(5):1466–1472.